

Dental CAM Library for **exocad**[®] Models

Installation Instructions

Tibase – Premill – Ang - Bars



Dental CAM
Digital Solutions

1. Download the folder *DentalCAM_library_mod_creator_ver3.1*
2. Unzip *DentalCAM_library_mod_creator_ver3.1*
3. Open the Folder *DentalCADApp/ library/modelcreator/implants* in your Exocad installation
4. Copy the folders you unzip
5. And paste into *DentalCADApp/ library/modelcreator/implants*



- 1- Descargue a su escritorio la carpeta *DentalCAM_library_mod_creator_ver3.1*
- 2- Unzip *DentalCAM_library_mod_creator_ver3.1*
- 3- Abra la carpeta *DentalCADApp/ library/modelcreator/implants* en su instalacion de Exocad
- 4- Copie las carpetas que descomprimio
- 5- Ahora pegue estas carpetas en *DentalCADApp/ library/modelcreator/implants*

